























ESAND THAILAND CODING & AI ACADEMY

โครงการวิจัยโมเดลระบบนิเวศการเรียนรู้ที่บูรณาการ CODING & AI **สำหรับเยาวชน** Model of Learning Ecosystem Platform integrate with Coding & Al for Youth



โครงการย่อยที่ 6

การพัฒนาเยาวชนเพื่อเข้าสู่วิชาชีพขั้นสูงด้าน Coding & Al ร่วมกับ Coding Entrepreneur & Partnership: Personal Al

BiTNet: AI for Ultrasound Image Classification

ผศ.ดร.ธนพงศ์ อินทระ ผู้เชี่ยวชาญด้าน Computer Vision



























E SAN THAILAND โครงการวิจัยโมเดลระบบนิเวศการเรียนรู้ที่บูรณาการ CODING & AI สำหรับเยาวชน Model of Learning Ecosystem Platform integrate with Coding & Al for Youth

Model Development




EfficientNet (Base Model)




BiTNet

CODING & AI ACADEMY Model of Learning Ecosystem Platform integrate with Coding & Al for Youth



EfficientNet (Base Model)



BiTNet

Biliary Tract Network



EfficientNet






Figure 2. Model Scaling. (a) is a baseline network example; (b)-(d) are conventional scaling that only increases one dimension of network width, depth, or resolution. (e) is our proposed compound scaling method that uniformly scales all three dimensions with a fixed ratio.



Performance Comparison of Base Models








Auto Pre-screening


2 Applications



1st Application



1st Application






Auto Pre-screening


#images identified as abnormal #submitted images

Comparison between workload reductionrate and false negative rate when variesthresholds of the model.








Auto Pre-screening




2nd Application








2nd Application





E SAU THAILAND โครงการวิจัยโมเดลระบบนิเวศการเรียนรู้ที่บูรณาการ co CODING & AI ACADEMY | Model of Learning Ecosystem Platform integra

Assisting tool







Group1 Group2

unassisted

Data distribution (150 test images)

	FP-A	FP-B	FP-C	FP-D	FP-E
AB01	1	1	1		
AB02	1	1	1		
AB03	1	1	1		
AB04	1	1	1	1	
AB05	1	1	1		
AB06	1	1	1		
AB07	1	1	1		
AB081	1	1	1		
AB082	1	1	1		
AB083	1	1	1		
AB09		2	1		
AB10			3		
AB11			1	2	
AB12				3	
Abnormal	11	12	14	6	0
Normal	22	24	28	12	21

Total: 150 images Abnormal: 43 images Normal: 107 images

Training session

Session 1:

Diagnose 150 test images

Washout Period of 4 Weeks


Session 2:


Diagnose 150 test images




5 general practitioners (GP's), 2 residence radiologists, 2 non-hepatobiliary radiologists and 2 hepatobiliary radiologists.





















Head








1,000 classes SoftMax

Freezed



Unfreezed

